

Product Datasheet

Human CCL4 protein (orb1216182)

Description	The Human CCL4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	CCL4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APMGSDPPTA CCFSYTARKL PRNFVVDYYE TSSLCSQPAV VFQTKRSKQV CADPSETWVQ EYVYDLELN (69)
Purity	98%
Source	Yeast
MW	7.8 kDa
Entrez	6351
Expiration Date	6 months from date of receipt.