

Product Datasheet

Canine CCL4 protein (orb1216183)

Description	The Canine CCL4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	CCL4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APMGSDPPTA CCFSYTLRKI PRNFVADYFE TSSLCSQPAV VFQTRRGRQV CANPSEPWWQ EYMDDLELN (69)
Purity	98%
Source	Yeast
MW	7.8 kDa
Entrez	448786
Expiration Date	6 months from date of receipt.