

## Product Datasheet

### Mouse B7-H2 protein (orb1216197)

<b>Description</b>	The Mouse B7-H2 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
<b>Target</b>	B7-H2
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	ETEVGAMVGS NVVLSCIDPH RRHFNLGLY VYWQIENPEV SVTYLPHYKS PGINVDSSYK NRGHLSLDSM KQGNFSLYLK NVTPQDTQEF TCRVFMNTAT ELVKILEEVV RLRVAANFST PVISTSDSSN PGQERTYTCM SKNGYPEPNL YWINTTDNSL IDTALQNNTV YLNKLGLYDV ISTLRLPWTS HGDVLCVEN VALHQNITSI SQAESFTGNN TKNPQETHNN ELKVLVPVLA VLAAAAFVSF IYRRTRPHR SYTGPKTVQL ELTDHA (276) [His185Arg]
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	31.0 kDa
<b>Entrez</b>	<a href="#">50723</a>
<b>Expiration Date</b>	6 months from date of receipt.