

Product Datasheet

Porcine BAFF protein (orb1216208)

| | |
|-------------------------|--|
| Description | The Swine BAFF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ... |
| Target | BAFF |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | AAQDAEETVT QDCLQLIADS DMPTIRKGAY TFVPWLLSFK RGRALEEKEN KIVVKETGYF FIYGQVLYTD NTFAMGHLIQ RKKVHVFGDE LSLVTLFRCI QNMPETLPNN SCYSAGIAKL EEGDELQLAI PREDAKISRD GDGTFFGALK LL (152) |
| Purity | 98% |
| Source | Yeast |
| MW | 17.1 kDa |
| Entrez | 100038026 |
| Expiration Date | 6 months from date of receipt. |