

## Product Datasheet

### Mouse IFN beta protein (orb1216213)

<b>Description</b>	The Mouse IFN beta yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. ...
<b>Target</b>	IFN beta
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	IN YKQLQLQE RTNIRKCQEL LEQLNGKINL TYRADFKIPM EMTEKMQKSY TAFAIQEMLQ NVFLVFRNNF SSTGWNETIV VRLDELHQQ TVFLKTVLEE KQEERLTWEM SSTALHLKSY YWRVQRYLKL MKYNSYAWMV VRAEIFRNFL IIRRLTRNFQ N (161)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	19.7 kDa
<b>Entrez</b>	<b>15977</b>
<b>Expiration Date</b>	6 months from date of receipt.