

## Product Datasheet

### Canine IL-1 alpha protein (orb1216214)

<b>Description</b>	The Canine IL-1 alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
<b>Target</b>	IL-1 alpha
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	SVAYNFHNNKEYNYIRIIKSQFILNDNLNQSIVRQTGGNYLMTAALQNLDDAVKF DMGAYTSEDSKLPVTLRISKTRLFVSAQNEDEPVLLKEMPETPKTIRDETNLFFW ERHGSKHYFKSVAQPKLFATQERKLVHMARGQPSITDFRLLLETQP (157)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	18.2 kDa
<b>Entrez</b>	<a href="#">403782</a>
<b>Expiration Date</b>	6 months from date of receipt.