

Product Datasheet

Canine CCL5 protein (orb1216224)

Description	The Canine CCL5 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	CCL5
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SPYASDTTPC CFAYISGRLP FTHVQEYFYT SSKCSMPAVV FVTRKHRQVC ANPQKKWVRE YINSLEMS (68)
Purity	98%
Source	Yeast
MW	7.9 kDa
Entrez	403522
Expiration Date	6 months from date of receipt.