

Product Datasheet

Mouse TNFSF11 protein (orb1216229)

Description	The Mouse TNFSF11 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	TNFSF11
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	FSGAPAMMEG SWLDVAQRGK PEAQPFAHLT INAASIPSGS HKVTLSSWYH DRGWAKISNM TLSNGKLRVN QDGFYYLYAN ICFRHHETSG SVPTDYLQLM VYVVKTSIKI PSSHNLKGG STKNWSGNSE FHFYSINVGG FFKLRAGEEI SIQVSNPSSL DPDQDATYFG AFKVQDID (178)
Purity	98%
Source	Yeast
MW	19.8 kDa
Entrez	21943
Expiration Date	6 months from date of receipt.