

Product Datasheet

Equine IL-1ra protein (orb1216234)

Description	The Equine IL-1 Receptor Antagonist yeast-derived recombinant protein is not tagged and is naturally...
Target	IL-1ra
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	HPLGKRPCCKMQAFRIWDVFNQKTFYMRNNQLVAGYLQESNTKLQEKIDVVPIEPD ALFLGLHGRKLCLACVKSGDEIRFQLEAVNITDLSKNKEENKRFTFIRSNSGPTTS FESAACPGWFLCTAQEADRPVSLTNKPKESFMVTKFYLQEDQ (152)
Purity	98%
Source	Yeast
MW	17.4 kDa
Entrez	100034236
Expiration Date	6 months from date of receipt.