

Product Datasheet

Mouse IL-6 protein (orb1216250)

Description	The Mouse IL-6 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	IL-6
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	FPTSQVRRGDFTE D T T P N R P V Y T T S Q V G G L I T H V L W E I V E M R K E L C N G N S D C M N N D D A L A E N N L K L P E I Q R N D G C Y Q T G Y N Q E I C L L K I S S G L L E Y H S Y L E Y M K N N L K D N K K D K A R V L Q R D T E T L I H I F N Q E V K D L H K I V L P T P I S N A L L T D K L E S Q K E W L R T K T I Q F I L K S L E E F L K V T L R S T R Q T (187)
Purity	98%
Source	Yeast
MW	21.7 kDa
Entrez	16193
Expiration Date	6 months from date of receipt.