

## Product Datasheet

### Mouse IL-1 beta protein (orb1216263)

<b>Description</b>	The Mouse IL-1 beta yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....
<b>Target</b>	IL-1 beta
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	VPIRQLHYRLRDEQQKSLVLSDPYELKALHLNGQNINQQVIFSMSFVQGEPNSDK IPVALGLKGKNLYLSCVMKDGTPTLQLESVDPKQYPKKKMEKRFVFNKIEVKSKV EFESAEFPNWIYSTSQAEHKPVFLGNNSGQDIIDFTMESVSS (152)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	17.4 kDa
<b>Entrez</b>	<b>16176</b>
<b>Expiration Date</b>	6 months from date of receipt.