

Product Datasheet

Mouse IL-17A protein (orb1216267)

Description	The Mouse IL-17A yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-17A
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	AAIIPQSSACPNTAKDFLQNVKVNKVFNSLGAKVSSRRPSDYLNIRSTSPWTLH RNEDPDRYPSVIWEAQCRHQRCVNAEGKLDHMHMNSVLIQQEILVLKREPESCPF TFRVEKMLVGVGCTCVASIVRQAA (133)
Purity	98%
Source	Yeast
MW	15.0 kDa
Entrez	16171
Expiration Date	6 months from date of receipt.