

Product Datasheet

Bovine IL-1ra protein (orb1216268)

Description	The Bovine IL-1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-1ra
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	HPLGKRRCEM QAFRIWDVNQ KIFYLRNNQL VAGYLQGPNT KLEEKIDVVP IEPHTMFLGI HGGKLCLACV KSGDEIKLKL EAVNITDLNQ NREQDKRFAF IRFDNGPTTS FESAACPGWF LCTSLEADQP VGLTNMPTEA LKVTKFYFQQ D (151)
Purity	98%
Source	Yeast
MW	17.2 kDa
Entrez	281860
Expiration Date	6 months from date of receipt.