

Product Datasheet

Guinea pig IL-8 protein (orb1216275)

| | |
|-------------------------|---|
| Description | The Guinea Pig IL-8 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free.... |
| Target | IL-8 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | MVVTKLVSEL RCQCIKIHTT PFHPKFIKEL KVIESGPRCA NSEIIVKLSL NRQLCLDPKK KVVQDVVSMF LKRTESQDS (79) |
| Purity | 98% |
| Source | Yeast |
| MW | 9.1 kDa |
| Entrez | 100379599 |
| Expiration Date | 6 months from date of receipt. |