

Product Datasheet

Bovine IL-16 protein (orb1216279)

Description	The Bovine IL-16 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-16
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	ATTDLNSSTDSSGSASVTS DVSIESAEATVCTVTLEKTSAGLGFSLEGGKGS LHG DKPLTVNRIFKGLASEQSDTVQPGDEIVHLAGTAMQGLTRFEAWNIIKALPDGPV TIVLRRKSLQSKGTPAAGDP (130)
Purity	98%
Source	Yeast
MW	13.3 kDa
Entrez	506314
Expiration Date	6 months from date of receipt.