

Product Datasheet

Mammal IL-8 protein (orb1216282)

Description	The Dolphin IL-8 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-8
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	AVLSRMTSEL RCQCINIHST PFHPKFIREL RVIESGPHCE NSEIIVKLVN GKEVCLNPKE KVVQKVVQIF LKRAEKKDP (79)
Purity	98%
Source	Yeast
MW	9.2 kDa
Entrez	not available
Expiration Date	6 months from date of receipt.