

## Product Datasheet

### Mammal TNF alpha protein (orb1216283)

<b>Description</b>	The Dolphin TNF alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
<b>Target</b>	TNF alpha
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	LRSSSKTSSN KPVAHVVANL STQGQLRWLN TYANTLLANS VKLEDNQLVV PTDGLYLIYS QVLFQGGCP STHLFLHTI SRIAVSYPTK VNLLSAIKSP CQRETPEGAE AKPWYEPIYQ GGVFQLEKGD RLSAEINLPD YLDAESGQV YFGIIL (157) [Ser89Thr]
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	17.4 kDa
<b>Entrez</b>	not available
<b>Expiration Date</b>	6 months from date of receipt.