

Product Datasheet

Bovine IL-18 protein (orb1216286)

Description	The Bovine IL-18 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-18
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	HFGKLEPKLS IIRNLNDQVL FINQGNQPVF EDMPDSDCSD NAPQTIFIY MYKDSLTRGL AVTISVQCKK MSTLSCENKI VSFKEMNPPD NIDNEESDII FFQRSVPGHD DKIQFESSLY EGYFLACKKE NDLFKLILKK QDDNRDKSVM FTVQNQN (157) [Lys121Glu]
Purity	98%
Source	Yeast
MW	18.1 kDa
Entrez	281249
Expiration Date	6 months from date of receipt.