

## Product Datasheet

### Porcine IL-17F protein (orb1216297)

<b>Description</b>	The Swine IL-17F yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
<b>Target</b>	IL-17F
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	ARKILKPGDP ALCPPLEDHS VRVDIRILSQ NQGVPLSHDF QNRSSSPWDY NITRDPHRFP SEIAEARCRH TGCINAQQQE DSSMNSVPIQ QEFLVLRREP QGCSRSFRLE KVRVTVGCTC VTPMVRYVRG LRAAGQLG (138)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	15.5 kDa
<b>Entrez</b>	<b>100037933</b>
<b>Expiration Date</b>	6 months from date of receipt.