

Product Datasheet

Canine IFN gamma protein (orb1216301)

Description	The Canine IFN gamma yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IFN gamma
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QAMFFKEIEN LKEYFNASNP DVSDGGSLFV DILKKWREES DKTIIQSQIV SFYLKLFDNF KDNQIIQRSM DTIKEDMLGK FLNSSTSKRE DFLKLIQIPV NDLQVQRKAI NELIKVMNDL SPRSNLRKRK RSQNLFRGRR ASK (143)
Purity	98%
Source	Yeast
MW	16.7 kDa
Entrez	403801
Expiration Date	6 months from date of receipt.