

Product Datasheet

Porcine IL-22 protein (orb1216304)

Description	The Swine IL-22 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-22
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	VPITHHCKLD QSNFQQPYIT NRTFTLAQEA SLADNNTDVR LIGNNLFQGV NMRERCYLVK QVLNFTLEEV LFPNSDRFHP YMQEVASFLD SLSKKLSQCR IKGDDQHIQR NVNFKDTVK KLGESGEIKV IGELYLLFMA LKNECTLPGH SWKMDN (156) [Gln51/52]
Purity	98%
Source	Yeast
MW	18.0 kDa
Entrez	595104
Expiration Date	6 months from date of receipt.