

Product Datasheet

Canine IL-21 protein (orb1216306)

| | |
|-------------------------|--|
| Description | The Canine IL-21 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th... |
| Target | IL-21 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | HKSSFQEQDL LLIRMRQLID IVDQLKNYVN DLDPEPLPAP EDVKRHCERS AFSCFQKVQL KAANTGGNEQ IINVLTKQLK RKLPTNAGR RQKHRPACPS CDSYEKAPPK EFLERLKSLI QKMIHQHLS (129) |
| Purity | 98% |
| Source | Yeast |
| MW | 14.9 kDa |
| Entrez | 442935 |
| Expiration Date | 6 months from date of receipt. |