

Product Datasheet

Rabbit IL-1 beta protein (orb1216314)

| | |
|-------------------------|--|
| Description | The Rabbit IL-1 beta yeast-derived recombinant protein is not tagged and is naturally endotoxin-free... |
| Target | IL-1 beta |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | AVRSLHCRLQ DAQQKSLVLS GTYELKALHL NAENLNQQVW FSMSFVQGE E SNDKIPVALG LRGKNLYLSC VMKDDKPTLQ LESVDPNRY P KKKMEKRFVF NKIEIKDKLE FESAQFPNWY ISTSQTEYMP VFLGNSSGGQ DLIDFSMEFV SS (152) |
| Purity | 98% |
| Source | Yeast |
| MW | 17.4 kDa |
| Entrez | 100008990 |
| Expiration Date | 6 months from date of receipt. |