

Product Datasheet

Porcine IL-6 protein (orb1216317)

| | |
|-------------------------|--|
| Description | The Swine IL-6 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ... |
| Target | IL-6 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | GRLEEDAKGDATSDKMLFTSPDKTEELIKYILGKISAMRKEMCEKYKENSKEVL AENNLNLPKMAEKDGCFSQGFNQETCLMRITTTGLVEFQIYLDYLQKEYESNKG VEAVQISTKALIQTLRQKGNPKATTPNPTTNAAGLLDKLQSQNEWMKNTKIILIL RSLEDFLQFSLRAIRIM (183) |
| Purity | 98% |
| Source | Yeast |
| MW | 20.9 kDa |
| Entrez | 399500 |
| Expiration Date | 6 months from date of receipt. |