

## Product Datasheet

### Equine TNF alpha protein (orb1216318)

|                         |  |
|-------------------------|--|
| <b>Description</b>      | The Equine TNF alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...  |
| <b>Target</b>           | TNF alpha  |
| <b>Form/Appearance</b>  | Lyophilized  |
| <b>Storage</b>          | -20°C  |
| <b>Note</b>             | For research use only  |
| <b>Protein Sequence</b> | LRSSSRTPSD KPVAVVAVNP QAEGQLQWLS<br>GRANALLANG VKLTDNQLVV PLDGLYLIYS<br>QVLFKQGQCP STHVLLHTI SRLAVSYPSK<br>VNLLSAIKSP CHTESPEQAE AKPWYEPIYL<br>GGVFQLEKGD QLSAEINQPN YLDFAESGQV<br>YLGIIAL (157) [Phe152Leu] |
| <b>Purity</b>           | 98%  |
| <b>Source</b>           | Yeast  |
| <b>MW</b>               | 17.1 kDa   |
| <b>Entrez</b>           | <b>100033834</b>   |
| <b>Expiration Date</b>  | 6 months from date of receipt.   |