

Product Datasheet

Rabbit IFN gamma protein (orb1216319)

Description	The Rabbit IFN gamma yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IFN gamma
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QDTLTRETEHLKAYLKANTSDVANGGPLFLNLRNWKEESDNKIIQSQIVSFYFKL FDNLKDHEVIKKSMESEKEDIFVKFFNSNLTKMDDFQNLTRISVDDRLVQRKAVS ELSNVLNFLSPKSNLKKRKRSQLFRGRRASKY (144)
Purity	98%
Source	Yeast
MW	16.9 kDa
Entrez	100008602
Expiration Date	6 months from date of receipt.