

Product Datasheet

Feline IFN gamma protein (orb1216320)

Description	The Feline IFN gamma yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IFN gamma
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QAMFFKEIEELKGYFNASNPDVADGGSLFVDILKNWKEESDKTIIQSQIVSFYMKM FENLKDDDDQRIQRSMDTIKEDMLDKLLNTSSSKRDDLKLIQIPVNDLQVQRKAI NELFKVMNDLSPRSNLRKRKRSQNLFRGRRASK (144)
Purity	98%
Source	Yeast
MW	17.0 kDa
Entrez	493965
Expiration Date	6 months from date of receipt.