

Product Datasheet

Porcine IFN gamma protein (orb1216323)

Description	The Swine IFN gamma yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....
Target	IFN gamma
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QAPFFKEITI LKDYFNASTS DVPNGGPLFL EILKNWKEES DKKIIQSQIV SFYKFFEIF KDNQAIQRSM DVIKQDMFQR FLNGSSGKLN DFEKLIKIPV DNLQIQRKAI SELIKVMNDL SPRSNLRKRK RSQTMFQGQR ASK
Purity	98%
Source	Yeast
MW	16.8 kDa
Entrez	396991
Expiration Date	6 months from date of receipt.