

Product Datasheet

Feline IL-1 beta protein (orb1216337)

Description	The Feline IL-1 beta yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IL-1 beta
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	AAIQSQDYTFRDISQKSLVLSGSYELRALHLNGQNMNQVVFRMSFVHGEENSK KIPVVLCIKKNNLYLSCVMKDGKPTLQLEMLDPKVYPKKKMEKRFVFNKTEIKGNV EFESSQFPNWIYSTSQAEEMPVFLGNTKGGQDITDFIMESAS
Purity	98%
Source	Yeast
MW	17.4 kDa
Entrez	768274
Expiration Date	6 months from date of receipt.