

## Product Datasheet

### Canine IL-1 beta protein (orb1216339)

|                         |   |
|-------------------------|---|
| <b>Description</b>      | The Canine IL-1 beta yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...   |
| <b>Target</b>           | IL-1 beta   |
| <b>Form/Appearance</b>  | Lyophilized   |
| <b>Storage</b>          | -20°C   |
| <b>Note</b>             | For research use only   |
| <b>Protein Sequence</b> | AAMQSV DCKLQDISHKYLVL SNSYELRALHLNGENVNKQVVFHMSFVHGDESN<br>NKIPVVLGIKQKNLYLSCVMKDGKPTLQLEKVDPKVYPKRKMEKRFVFNKIEIKNT<br>VEFESSQYPNWIYISTSQVEGMPVFLGNTRGGQDITDFTMEFSS |
| <b>Purity</b>           | 98%   |
| <b>Source</b>           | Yeast   |
| <b>MW</b>               | 17.5 kDa  |
| <b>Entrez</b>           | <a href="#">403974</a>  |
| <b>Expiration Date</b>  | 6 months from date of receipt.  |