

Product Datasheet

Porcine IL-15 protein (orb1216364)

Description	The Swine IL-15 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-15
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	TWQHVISDLK KIEDLIRSIH MDATLYTESD AHPNCKVTAM KCFLLELRVI LQESRNSDIS DTVENLILA NSSLSSIEYK TESGCKECEE LEEKNINEFL KSFHIVQMF INPS
Purity	98%
Source	Yeast
MW	13.1 kDa
Entrez	397683
Expiration Date	6 months from date of receipt.