

## Product Datasheet

### Bovine IL-6 protein (orb1216379)

<b>Description</b>	The Bovine IL-6 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
<b>Target</b>	IL-6
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	GPLGEDFKND TTPGRLLLT PEKTEALIKR MVDKISAMRK EICEKNDECE SSKETLAENK LNLPKMEEKD GCFQSGFNQA ICLIRTTAGL LEYQIYLDYL QNEYEGNQEN VRDLRKNIRT LIQILKQKIA DLITTPATNT DLLEKMQSSN EWWKNAKIIIL ILRNLENFLQ FSLRAIRMK
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	20.7 kDa
<b>Entrez</b>	<b>280826</b>
<b>Expiration Date</b>	6 months from date of receipt.