

Product Datasheet

Bovine IL-15 protein (orb1216384)

Description	The Bovine IL-15 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-15
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	NWQYVINDLKTIEHLIQSIHMDATLYTESDAHPNCKVTAMQCFLLELRVILHESKN ATIYEIENLTMLANSNLSSIENKTELGCKECEEELEEKSIKEFLKSFVHIVQMFINTS
Purity	98%
Source	Yeast
MW	13.2 kDa
Entrez	281248
Expiration Date	6 months from date of receipt.