

Product Datasheet

Equine IL-4 protein (orb1216390)

Description	The Equine IL-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GCKYDIT LQEIIKTLNN LTDGKGKNSC MELTVADAFAGPKNTDGKEICRAAKVLQQL YKRHDRSLIKECLSGLDNRNLKGMANGTCCT VNEAKKSTLKDFLERLKTIMKEKYSKC
Purity	98%
Source	Yeast
MW	12.8 kDa
Entrez	100034225
Expiration Date	6 months from date of receipt.