

Product Datasheet

CCL3 antibody (Biotin) (orb1216670)

Description	The Zebrafish CCL3 polyclonal antibody is labeled and has been qualified for use in ELISA. It is the...
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Zebrafish CCL3
Target	CCL3
Storage	4 degrees C.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.