

Product Datasheet

CXCL10 antibody (Biotin) (orb1216700)

Description	The Ferret/Mink CXCL10 (IP-10) polyclonal antibody is labeled and has been qualified for use in ELIS...
Species/Host	Rabbit
Reactivity	Ferret
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Ferret CXCL10
Target	CXCL10
Storage	4 degrees C.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.