

Product Datasheet

CCL4 antibody (Biotin) (orb1216718)

Description	The Cynomolgus Monkey CCL4 (MIP-1 beta) polyclonal antibody is labeled and has been qualified for us...
Species/Host	Rabbit
Reactivity	Primate
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Cynomolgus Monkey CCL4
Target	CCL4
Storage	4 degrees C.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.