

Product Datasheet

CCL11 antibody (Biotin) (orb1216736)

Description	The Human CCL11 (Eotaxin-1) polyclonal antibody is labeled and has been qualified for use in ELISA. ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human CCL11
Target	CCL11
Storage	4 degrees C.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.