

## Product Datasheet

### CXCL3 antibody (Biotin) (orb1216815)

<b>Description</b>	The Rabbit CXCL3 (GRO gamma) polyclonal antibody is labeled and has been qualified for use in ELISA....
<b>Species/Host</b>	Gallus
<b>Reactivity</b>	Rabbit
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Recombinant Rabbit CXCL3
<b>Target</b>	CXCL3
<b>Storage</b>	4 degrees C.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Expiration Date</b>	12 months from date of receipt.