

Product Datasheet

CCL21 antibody (Biotin) (orb1216911)

Description	The Rabbit CCL21 (6Ckine) polyclonal antibody is labeled and has been qualified for use in ELISA. It...
Species/Host	Gallus
Reactivity	Rabbit
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Rabbit CCL21
Target	CCL21
Storage	4 degrees C.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.