

## Product Datasheet

### CCL21 antibody (orb1216912)

<b>Description</b>	The Rabbit CCL21 (6Ckine) polyclonal antibody is unlabeled and has been qualified for use in ELISA a...
<b>Species/Host</b>	Gallus
<b>Reactivity</b>	Rabbit
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	Recombinant Rabbit CCL21
<b>Target</b>	CCL21
<b>Storage</b>	4 degrees C.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Expiration Date</b>	12 months from date of receipt.