

## Product Datasheet

### ARV-110 (orb1218412)

|                        |                                                       |
|------------------------|-------------------------------------------------------|
| <b>Description</b>     | ARV-110                                               |
| <b>Storage</b>         | Storage temperature: -20°C. Stability: $\geq 2$ years |
| <b>Note</b>            | For research use only                                 |
| <b>Purity</b>          | > 98% (HPLC)                                          |
| <b>CAS Number</b>      | 2222112-77-6                                          |
| <b>Expiration Date</b> | 12 months from date of receipt.                       |