

## Product Datasheet

### Beta-boswellic acid (orb1219414)

|                        |   |
|------------------------|---|
| <b>Description</b>     | Beta-boswellic acid                                 |
| <b>Storage</b>         | Storage temperature: -20°C. Stability:<br>≥ 2 years |
| <b>Note</b>            | For research use only                               |
| <b>Purity</b>          | > 98% (HPLC)  |
| <b>CAS Number</b>      | 631-69-6  |
| <b>Expiration Date</b> | 12 months from date of receipt.                     |