

## Product Datasheet

### ASLAN003 (orb1219816)

|                        |   |
|------------------------|---|
| <b>Description</b>     | ASLAN003  |
| <b>Storage</b>         | Storage temperature: -20°C. Stability: $\geq$ 2 years |
| <b>Note</b>            | For research use only                                 |
| <b>Purity</b>          | > 98% (HPLC)  |
| <b>CAS Number</b>      | 1035688-66-4  |
| <b>Expiration Date</b> | 12 months from date of receipt.                       |