

## Product Datasheet

### D-Lactose monohydrate (orb1221238)

|                        |                                                     |
|------------------------|-----------------------------------------------------|
| <b>Description</b>     | D-Lactose monohydrate                               |
| <b>Storage</b>         | Storage temperature: -20°C.<br>Stability: ≥ 2 years |
| <b>Note</b>            | For research use only                               |
| <b>Purity</b>          | > 98% (HPLC)                                        |
| <b>CAS Number</b>      | 64044-51-5                                          |
| <b>Expiration Date</b> | 12 months from date of receipt.                     |