

## Product Datasheet

### alpha-Arbutin (orb1221268)

|                        |                                                  |
|------------------------|--------------------------------------------------|
| <b>Description</b>     | alpha-Arbutin                                    |
| <b>Storage</b>         | Storage temperature: -20°C. Stability: ≥ 2 years |
| <b>Note</b>            | For research use only                            |
| <b>Purity</b>          | > 98% (HPLC)                                     |
| <b>CAS Number</b>      | 84380-01-8                                       |
| <b>Expiration Date</b> | 12 months from date of receipt.                  |