



Product Datasheet

Sp100 antibody (orb122509)

Description	Polyclonal Antibody pABahSp100 (209-227) # 2 to human Sp100.
Species/Host	Rabbit
Conjugation	Unconjugated
Immunogen	The antibody binds to a peptide corresponding to amino acids 209-227 (CETEQINAKRKDTTSDKDD) of the human Sp100 protein.
Form/Appearance	liquid (Supplied as serum)
Storage	store at -20 °C
Note	For research use only
Clonality	Recombinant
Dilution Range	Indirect Immunofluorescence, concentration: 1:1000 to 1:2000. Western Blot, concentration: 1:500.
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of using Sp1...



Immunofluorescence analysis of using Sp1...