



Product Datasheet

PML antibody (orb122542)

Description	Polyclonal Antibody pABahPML (460-479) # 3 to human PML.
Species/Host	Rabbit
Conjugation	Unconjugated
Immunogen	The antibody binds to a peptide corresponding to amino acids 460-479 (KGPSYGEDVSNTTTQKRKC) of the human PML protein. No cross reaction with mouse PML protein was observed in Western blots. The immunogen peptide sequence is not conserved in rat, chimpanzee, dog, or bovine PML proteins.
Form/Appearance	liquid (Supplied as serum)
Storage	store at -20 °C
Note	For research use only
Clonality	Recombinant
Dilution Range	Indirect Immunofluorescence 1:250. Western Blot 1:500
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of using PML...



Western blot analysis of using PML antib...