



Product Datasheet

TNNT2 antibody (orb12268)

Description	Goat polyclonal antibody to TNNT2
Species/Host	Goat
Reactivity	Bovine, Canine, Feline, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	cardiac troponin T
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 40kDa band observed in Human, fetal Mouse and adult Rat and Pig Heart lysates (calculated MW of 35.6kDa according to Human NP_000355.2 and 35.7kDa according to Rat NP_036808.1). The observed molecular weight corresponds to earlier findings in liter
Protein Sequence	RNRINDNQKVSKT
Clonality	Polyclonal
MW	35.6; 34.6; 34.3; 34.0
Entrez	7139
Dilution Range	ELISA: 1:128000, WB: 0.03-0.1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of staining of Rat...



Western blot analysis of staining of Pig...