

Product Datasheet

PDL1 ELISA Matched Antibody Pair (Risk Free) (orb1227433)

Description	PDL1 ELISA Matched Antibody Pair (Risk Free)
Species/Host	Mouse and Llama
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	PD-L1 antibodies were raised against the extracellular domain of human PD-L1.
Preservatives	PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	Antibody 1 mg/mL
Storage	Mouse Monoclonal: Antibodies can be stored at 4°C for three months and -20°C, stable for up to one year. Single Domain: Antibody should be stored in working aliquots at -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Isotype	IgG1 and sdab
Clonality	Monoclonal
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)